



## ATW0080

### APTAMER INFORMATION

### Brain Natriuretic Peptide

#### Description:

- **Identifiers:** BNP D2A4F
- **Number of DNA nucleotides:** 32 bases
- **Target for selection:** Brain natriuretic peptide, BNP (BACHEM, Cat # H-9060.0500)


Aptamer was selected at Base Pair Biotechnologies from a randomized 32-mer DNA library against the target protein. Proprietary methods were then used to select the aptamer.

#### Aptamer folding instruction before use:

Once the aptamer is in its working concentration, it needs to be heated to 85-95 °C for 5 minutes, and then cooled to room temperature (~ 15 minutes) before use.

#### Validation data by Microscale Thermophoresis (MST) analysis:

- **Buffer used for validation:** 1X PBS, pH = 7.4, 1 mM MgCl<sub>2</sub>, 0.05% Tween
- **Average  $K_d$ :** 60.4 ± 10.9 nM



**Figure 1.** Aptamer-BNP binding. The fraction of bound aptamer is plotted versus the titrated target concentration,  $R^2 = 0.986$ .